


SIMPLY EASY




Simply see test results.

Large, easy-to-read display makes it simple to see test results.



Simply use the entire end.

Place a small drop of blood anywhere along the end instead of a tiny spot.



Simply take just one strip.

Spill-resistant SmartPack® vial lets you take one and spill none.



Simply lance with less pain.

The Accu-Chek® Softclix lancing device provides personalized comfort and unique technology to minimize pain.



Simply more accurate.

Accu-Chek Guide test strips deliver advanced 10/10 accuracy for reliable results you can trust.¹

Additional information can be found at hcp.accu-chek.com

¹ 95% of the measured glucose results shall fall within 10 mg/dL of lab reference value for blood glucose concentrations <100 mg/dL and within 10% for blood glucose concentrations ≥100 mg/dL. Data on file.