

Edit an existing patient details

Step One: Patient Profile

From the Patient Dashboard, click on the Patient Settings drop down arrow and select Edit Patient.

The screenshot shows the Patient Dashboard for 'David Demo'. The patient's details are: Patient Number 12345, Date of Birth 02/01/2009, Diabetes Type TYPE 2, and Center NAM. The dashboard includes a 'Patient Settings' dropdown menu with options: Patient Profile, Edit Patient, Time Blocks, Glucose Target Ranges, Devices, and Deactivate Patient. The 'Edit Patient' option is highlighted with a dashed red box. Below the menu, the 'GLUCOSE LEVEL DISTRIBUTION' chart shows a target range of 70-180 mg/dL with 93% of readings within range. The 'STATUS' section shows 'Hypo Risk: High' and 'Mean Blood Glucose: Within target range (target range: 70-180 mg/dL)'. The 'VARIABILITY' section shows 'Variability: Low' at 14.6%.

Step Two: Edit Date

A new screen will appear, allowing you to edit patient details as necessary in Profile Type, Patient Info and Health Info tabs.

Click on **Save**.

The screenshot shows the 'Edit Patient Info' form. The form is divided into three tabs: PROFILE TYPE, PATIENT INFO, and HEALTH INFO. The 'PATIENT INFO' tab is selected. The form contains the following fields:

- First Name: David
- Last Name(s): Demo
- Date of Birth Optional: Jan, 02, 2009
- Healthcare ID: 12345